ACCURIS TM
SmartCheck™ DNA Ladders
□ PR4005 50 bp ladder □ Sample 1 x 25ūl □ PR4010 100 bp ladder □ 100 lanes 1 x 500ul □ PR4100 1kb ladder □ 500 lanes 5 x 500ul
Ready-to-use with loading dye added, 0.1mg/ml concentration Store at -20°C upon receipt.

SmartCheckTM DNA Ladders PR4005 50 bp ladder Sample PR4010 100 lanes 100 bp ladder 1 x 500ul PR4100 500 lanes 1kb ladder 5 x 500ul Ready-to-use with loading dye added, 0.1mg/ml concentration Store at -20°C upon receipt. PH: 908-769-5555 EM: info@accuris-usa.com

Accuris™ SmartCheck™ DNA Ladders

Description

Accuris SmartCheck DNA Ladders are produced from proprietary plasmids, digested to completion and purified in a multi-step chromatography process. Ladders are supplied in a ready-to-use format that includes loading buffer and tracking dyes.

- -Three sizes: 50bp, 100bp and 1Kb
- -Clear, distinct bands
- -Reference bands with increased intensity for easy size determination
- -Convenient, simply load 5µl per lane

PR4010 / PR4100 ladders are supplied in a buffer consisting of 10mM Tris-HCI (pH 8.0), 5mM EDTA, 12.5% glycerol, 0.008% bromophenol blue, 0.008% xylene cyanol.

PR4005 ladders are supplied in a buffer consisting of 10mM Tris-HCl (pH 8.0), 10mM Na-EDTA, 0.025% Orange G, 10% glycerol.

Storage

This product is shipped at ambient temperature. Upon receipt, ladders can be stored for up to 24 months at -20°C or at 4°C for up to 3 months.

Quality Control

The ladders are purified by a multistep chromatography method. Band sharpness is verified by agarose electrophoresis and DNA concentration is verified spectrophotometrically.

General Guidelines

SmartCheck DNA ladders are supplied ready-to-use. Gently pipette up and down to mix the ladders and add 5μ l to the desired well in the gel. If the ladders are frozen, allow time to thaw and mix thoroughly.

The masses of each DNA band per $5\mu l$ load are given on the back side of this insert.

Technical Support

For trouble-shooting and technical support, contact us by phone at 908 769-5555 or email info@accuris-usa.com.

Limitations of Use

For research purposes only. Not intended for therapeutic or diagnostic use. Accuris is not responsible for consequential or incidental damages, whether direct or indirect, resulting from use of this product.

Package contents and reordering

See front panel for package contents and reordering information.

See reverse side for images of each ladder with band sizes and masses. Images and be cut and placed on lab bench for quick reference.

Accuris offers a full line of PCR enzymes and master mixes. Visit www.accuris-usa.com for details.

Accuris™ SmartCheck™ DNA Ladders

Description

Accuris SmartCheck DNA Ladders are produced from proprietary plasmids, digested to completion and purified in a multi-step chromatography process. Ladders are supplied in a ready-to-use format that includes loading buffer and tracking dyes.

- -Three sizes: 50bp, 100bp and 1Kb
- -Clear, distinct bands
- -Reference bands with increased intensity for easy size determination
- -Convenient, simply load 5µl per lane

The ladders are supplied at a concentration of 0.1mg/ml, in a buffer consisting of 10mM Tris-HCl (pH 8.0), 5mM EDTA, 12.5% glycerol, 0.008% bromophenol blue, 0.008% xylene cyanol.

PR4005 ladders are supplied in a buffer consisting of 10mM Tris-HCl (pH 8.0), 10mM Na-EDTA, 0.025% Orange G, 10% glycerol.

Storage

This product is shipped at ambient temperature. Upon receipt, ladders can be stored for up to 24 months at -20°C or at 4°C for up to 3 months.

Quality Control

The ladders are purified by a multistep chromatography method. Band sharpness is verified by agarose electrophoresis and DNA concentration is verified spectrophotometrically.

General Guidelines

SmartCheck DNA ladders are supplied ready-to-use. Gently pipette up and down to mix the ladders and add 5μ l to the desired well in the gel. If the ladders are frozen, allow time to thaw and mix thoroughly.

The masses of each DNA band per $5\mu l$ load are given on the back side of this insert.

Technical Support

For trouble-shooting and technical support, contact us by phone at 908 769-5555 or email info@accuris-usa.com.

Limitations of Use

For research purposes only. Not intended for therapeutic or diagnostic use. Accuris is not responsible for consequential or incidental damages, whether direct or indirect, resulting from use of this product.


Package contents and reordering

See front panel for package contents and reordering information.

See reverse side for images of each ladder with band sizes and masses. Images and be cut and placed on lab bench for quick reference.

Accuris offers a full line of PCR enzymes and master mixes. Visit www.accuris-usa.com for details.


50 bp ladder - PR4005



2% Agarose Gel, 0.5x TAE Buffer


Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

100 bp ladder - PR4010



Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

1 kb ladder - PR4100



Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

Electrophoresis Products and Accessories



302nm UV wave length,

High quality, enhanced





MYVOLT "MINI


supply with color touch screen



100/200V output power supply


For additional product details, please visit www.accuris-usa.com or contact your local Benchmark Scientific dealer

50 bp ladder - PR4005




2% Agarose Gel, 0.5x TAE Buffer Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

100 bp ladder - PR4010



Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

1 kb ladder - PR4100



Accuris SmartCheck™ DNA Ladders www.accuris-usa.com

ACCURIS Electrophoresis Products and Accessories

Benchmark's Accuris Division offers a wide range of products to make electrophoresis easier, safer and smarter. Contact your local distributor for details to place an order.



fluoresces under UV or blue light



MY GEL



with a smart phone or tablet



contact your local Benchmark Scientific dealer.