




## SmartGlow™ Pre Stain, for Nucleic Acid Electrophoresis Gels, Concentration 20,000X

### E4500-PS

SmartGlow™ Pre Stain can be used the same way as Ethidium Bromide for agarose and polyacrylamide gel electrophoresis. The stain emits green fluorescence when bound to dsDNA, ssDNA, and RNA. SmartGlow™ stain has two excitation peaks at approximately 290nm (UV light) and 475nm (Blue light).



SmartGlow™ is for research use only


Ordering Information:  
E4500-PS 1.0m SmartGlow™ Pre Stain, for Nucleic Acid Electrophoresis Gels, 20,000X  
Storage: 4°C for 2 years



## SmartGlow™ Pre Stain, for Nucleic Acid Electrophoresis Gels, Concentration 20,000X

### E4500-PS

SmartGlow™ Pre Stain can be used the same way as Ethidium Bromide for agarose and polyacrylamide gel electrophoresis. The stain emits green fluorescence when bound to dsDNA, ssDNA, and RNA. SmartGlow™ stain has two excitation peaks at approximately 290nm (UV light) and 475nm (Blue light).



SmartGlow™ is for research use only

Ordering Information:  
E4500-PS 1.0m SmartGlow™ Pre Stain, for Nucleic Acid Electrophoresis Gels, 20,000X  
Storage: 4°C for 2 years



**Benchmark**  
*Scientific*  
P.O. Box 709 Edison, NJ 08818 USA  
Ph: 908-769-5555 EM: info@accuris-usa.com

SmartGlow™ Pre Stain is considered safer than Ethidium Bromide, and is non-carcinogenic as determined by the Ames-test, with negative results in mouse primary spermatocyte chromosomal aberration tests and mouse marrow chromophilous erythrocyte micronucleus tests. SmartGlow™ Pre Stain can be disposed as non-hazardous waste, and has passed the CCR Title 22 Hazardous Waste Test.

- Protocol:
- Prepare 100ml of agarose (heat as required) or 100ml polyacrylamide solution.
  - For agarose, allow the solution to cool to about 60°C to 70°C.
  - Add 5μl\* of SmartGlow™ stain to the gel solution.
  - Mix gently, make sure there are no air bubbles, pour the gel(s).
  - For best results, add 5μl of SmartGlow™ stain to each 100ml of running buffer.
  - Adding the stain to the running buffer increases sensitivity.
  - After the separation is complete, view the gel using a UV, blue light illuminator, or SmartDoc system.
- \* If DNA/background fluorescence is too high, use 2.5μl/100ml.

## Also available from Benchmark Scientific and Accuris Instruments:



### AGAROSE LE

- A1700 Benchmark Agarose LE, 25g
- A1701 Benchmark Agarose LE, 100g
- A1705 Benchmark Agarose LE, 500g



Includes two casting sets

- E1101 MyGel Mini Electrophoresis System, complete with gel casting sets

### SMARTBLUE™



- E4000 SmartBlue, blue light transilluminator



E5001-SDB SmartDoc Gel Imaging System with blue light illumination

E5001-SD SmartDoc Gel Imaging Enclosure (for use on existing transilluminators)



Contact us for more info or a price quote: [www.accuris-usa.com](http://www.accuris-usa.com) - [info@accuris-usa.com](mailto:info@accuris-usa.com) - 908-769-5555

## Also available from Benchmark Scientific and Accuris Instruments:



### AGAROSE LE

- A1700 Benchmark Agarose LE, 25g
- A1701 Benchmark Agarose LE, 100g
- A1705 Benchmark Agarose LE, 500g



Includes two casting sets

- E1101 MyGel Mini Electrophoresis System, complete with gel casting sets

### SMARTBLUE™



- E4000 SmartBlue, blue light transilluminator



E5001-SDB SmartDoc Gel Imaging System with blue light illumination

E5001-SD SmartDoc Gel Imaging Enclosure (for use on existing transilluminators)



Contact us for more info or a price quote: [www.accuris-usa.com](http://www.accuris-usa.com) - [info@accuris-usa.com](mailto:info@accuris-usa.com) - 908-769-5555